

SILICONE WOUND DRAINAGE SYSTEM

Silicone for great biocompatibility
and fast recovery



GENERAL
SURGICAL
SERIES

SILICONE WOUND DRAINAGE SYSTEM

Post-operative wound drainage can be safe and easy.



Made of medical grade silicone, FortuneMed® provides a complete line of C.W.V. (closed wound vacuum) Reservoirs, Drains, Trocars and Collection Bags for different surgical applications.

Why Silicone?

The advantages of silicone provide superior biocompatibility, flexibility and great comfort to patient.

For this reason, Fortune® silicone wound drain systems deliver extremely low tissue responses to wound, which minimize possible tissue-adhesion while healing and tissue-trauma while drain removal compare to drains of other elastomers, e.g. PVC and Latex.

- A Silicone C.W.V. Reservoir**
- B Collection Bag**
- C Trocar**
- D Silicone Flat Drain**
- E Round Drain**
- F Silicone Vac Drain**
- G Silicone Endo Drain**
- H Endo Connector**




*Data based on 102 female patients operated, 52 using PVC wound drainage system and 50 using Fortune® Silicone Wound Drain System. Drain pull-out time for measuring fluid drawn volume is 6.18 days for PVC system and 5.18 days for silicone system.

SILICONE C.W.V. RESERVOIR

Fortune® Silicone C.W.V. Reservoir Features:

- Medical grade silicone for superior biocompatibility.
- Low pressure design provides continuous gentle suction and minimizes tissue trauma.
- Volume marking provides easy measurement.
- Flat bottom design for easy positioning.
- Non-reflex valve.
- Hanger and bedside clip for secured attachment.
- Available with 4 sizes, 2 outlet directions.
- Optional extra inlet is available on 200 & 400 cc models to connect second tube.
- Bottom outlet models are capable to connect with Collection Bag.
- Sterilized double package.


A Silicone C.W.V. Reservoir



SPECIFICATIONS

C.W.V. Reservoir:

REF. No.	Size	Inlet	Description
2013-0610	100 cc	1	
2013-0615	150 cc		Top outlet, Bedside clip
2013-0620	200 cc	1+1	
2013-0640	400 cc		
2013-1640	400 cc	1+1	Top outlet, female luer w/cap Bedside clip
2013-0910	100 cc	1	
2013-0915	150 cc		Bottom outlet w/ large bore male connector. Pinch clamp & cap. Bedside clip
2013-0920	200 cc	1+1	
2013-0940	400 cc		



COLLECTION BAG

Fortune® Collection Bag Features:

- Large bore female connector for connecting with bottom outlet Reservoir.
- Non-reflex Valve.
- Volume marking provides easy measurement.
- Attached bag ties for secured bedside attachment.
- Increased protection compare to direct handling and disposal.
- Facilitates collected fluid sampling.
- Reduces bio-hazard for medical personnel.


A Collection Bag




SPECIFICATIONS

Collection Bag:


REF. No.	Size	Description
2013-9920	200 cc	- Large bore female connector
2013-9940	400 cc	- Pinch clamp & cap
2013-9980	800 cc	- Bag ties



2013-9920




2013-9940



2013-9980

Silicone C.W.V. Reservoir (Bottom outlet) + Collection Bag

REF. No.	Description
0910-9920	100 cc Reservoir, 200 cc Bag
0910-9940	100 cc Reservoir, 400 cc Bag
0915-9920	150 cc Reservoir, 200 cc Bag
0915-9940	150 cc Reservoir, 400 cc Bag
0920-9940	200 cc Reservoir, 400 cc Bag
0920-9980	200 cc Reservoir, 800 cc Bag
0940-9940	400 cc Reservoir, 400 cc Bag
0940-9980	400 cc Reservoir, 800 cc Bag



0910-9920

WOUND DRAINAGE TROCAR

Fortune® Wound Drainage Trocar features:

- Made of stainless steel.
- Sharp, smooth and burr free.
- Protection cover.
- Customized specification is available upon request.

A Trocar

SPECIFICATIONS

Trocar:

REF. No.	Size	Description
2015-9008	8 FR	
2015-9010	10 FR	
2015-9016	16 FR	Stainless steel, with Protection cover
2015-9019	19 FR	
2015-9020	20 FR	
2015-9021	21 FR	



2015-9021

SILICONE FLAT DRAIN

Fortune® Silicone Flat Drain Features:

- Medical grade silicone for superior biocompatibility.
- X-ray opaque effectiveness.
- Soft yet kink-resistant structure to minimize clog or obstruct.
- Clear connecting tube provides inspection and drainage visibility.
- Available with 4 sizes, 2 perforations types.
- Sterilized double package.
- Available with pre-connected Trocar: Silicone Flat Drain Kit.

- A Silicone Flat Drain**
- B Silicone Flat Drain Kit**
- C Silicone Flat Drain Set***
- D Silicone Flat Drain System***
- E Silicone Flat Drain Bag System***



SPECIFICATIONS

Flat Drain:

REF. No.	Size	Description
2014-0051	2 x 4.5 mm	
2014-0071	3 x 7 mm	3/4 perforation on 200 mm drain X-ray opaque effective
2014-0101	4 x 10 mm	
2014-0131	5 x 13 mm	
2014-0052	2 x 4.5 mm	
2014-0072	3 x 7 mm	Full perforation on 200 mm drain X-ray opaque effective
2014-0102	4 x 10 mm	
2014-0132	5 x 13 mm	




2014-OXX1 3/4 Perforation




2014-OXX2 Full Perforation

Flat Drain Kit:

REF. No.	Size	Description
2014-9051	2 x 4.5 mm	4.5 mm 3/4 Flat Drain, attached Trocar
2014-9071	3 x 7 mm	7 mm 3/4 Flat Drain, attached Trocar
2014-9101	4 x 10 mm	10 mm 3/4 Flat Drain, attached Trocar
2014-9131	5 x 13 mm	13 mm 3/4 Flat Drain, attached Trocar
2014-9052	2 x 4.5 mm	4.5 mm full Flat Drain, attached Trocar
2014-9072	3 x 7 mm	7 mm full Flat Drain, attached Trocar
2014-9102	4 x 10 mm	10 mm full Flat Drain, attached Trocar
2014-9132	5 x 13 mm	13 mm full Flat Drain, attached Trocar



2014-9XX1 3/4 Flat Drain with Trocar




2014-9XX2 Full Flat Drain with Trocar

SILICONE FLAT DRAIN

Flat Drain Set*:


REF. No.	Size	Description
0610-0051	2 x 4.5 mm	4.5 mm 3/4 Flat Drain, 100 cc Reservoir
0610-0071	3 x 7 mm	7 mm 3/4 Flat Drain, 100 cc Reservoir
0610-0101	4 x 10 mm	10 mm 3/4 Flat Drain, 100 cc Reservoir
0610-0131	5 x 13 mm	13 mm 3/4 Flat Drain, 100 cc Reservoir
0610-0052	2 x 4.5 mm	4.5 mm full Flat Drain, 100 cc Reservoir
0610-0072	3 x 7 mm	7 mm full Flat Drain, 100 cc Reservoir
0610-0102	4 x 10 mm	10 mm full Flat Drain, 100 cc Reservoir
0610-0132	5 x 13 mm	13 mm full Flat Drain, 100 cc Reservoir




0610-OXX1 3/4 Perforation Flat Drain Set

Flat Drain System *:


REF. No.	Size	Description
0610-9051	2 x 4.5 mm	4.5 mm 3/4 Flat Drain, Trocar, 100 cc Reservoir
0610-9071	3 x 7 mm	7 mm 3/4 Flat Drain, Trocar, 100 cc Reservoir
0610-9101	4 x 10 mm	10 mm 3/4 Flat Drain, Trocar, 100 cc Reservoir
0610-9131	5 x 13 mm	13 mm 3/4 Flat Drain, Trocar, 100 cc Reservoir
0610-9052	2 x 4.5 mm	4.5 mm full Flat Drain, Trocar, 100 cc Reservoir
0610-9072	3 x 7 mm	7 mm full Flat Drain, Trocar, 100 cc Reservoir
0610-9102	4 x 10 mm	10 mm full Flat Drain, Trocar, 100 cc Reservoir
0610-9132	5 x 13 mm	13 mm full Flat Drain, Trocar, 100 cc Reservoir



0610-9XX1 3/4 Flat Drain System



0610-OXX2 Full Perforation Flat Drain Set



0610-9XX2 Full Flat Drain System


SILICONE FLAT DRAIN

Flat Drain Bag System:


REF. No.	Size	Description
0910-905120	2 x 4.5 mm	4.5 mm 3/4 Flat Drain, Trocar, 100 cc Reservoir, 200 cc Bag
0910-907120	3 x 7 mm	7 mm 3/4 Flat Drain, Trocar, 100 cc Reservoir, 200 cc Bag
0910-910120	4 x 10 mm	10 mm 3/4 Flat Drain, Trocar, 100 cc Reservoir, 200 cc Bag
0910-913120	5 x 13 mm	13 mm 3/4 Flat Drain, Trocar, 100 cc Reservoir, 200 cc Bag
0910-905220	2 x 4.5 mm	4.5 mm full Flat Drain, Trocar, 100 cc Reservoir, 200 cc Bag
0910-907220	3 x 7 mm	7 mm full Flat Drain, Trocar, 100 cc Reservoir, 200 cc Bag
0910-910220	4 x 10 mm	10 mm full Flat Drain, Trocar, 100 cc Reservoir, 200 cc Bag
0910-913220	5 x 13 mm	13 mm full Flat Drain, Trocar, 100 cc Reservoir, 200 cc Bag

* Flat Drain Set, Flat Drain System and Flat Drain Bag System are also available with 150 cc, 200 cc or 400 cc Reservoir.

* Flat Drain Bag System is also available with 400 cc or 800 cc Collection Bag.



0910-9X0120 3/4 Flat Drain Bag System




0910-9X0220 Full Flat Drain Bag System

SILICONE ROUND DRAIN

Fortune® Silicone Round Drain Features:

- Medical grade silicone for superior biocompatibility.
- Clear tubing for visual drainage inspection.
- X-ray opaque line.
- Rigid and flexible structure for kinking resistance.
- Clear tubing provides inspection and drainage visibility.
- Available with 5 sizes, 2 perforations types.
- Sterilized double package.
- Available with pre-connected Trocar: Round Drain Kit.

- A Silicone Round Drain**
- B Silicone Round Drain Kit**
- C Silicone Round Drain Set***
- D Silicone Round Drain System***
- E Silicone Round Drain Bag System***



SPECIFICATIONS

Round Drain:

REF. No.	Size	Description
2015-0208	8 FR	
2015-0210	10 FR	
2015-0212	12 FR	- End perforation
2015-0214	14 FR	- 900 mm,
2015-0216	16 FR	- X-ray opaque line
2015-0218	18 FR	
2015-0220	20 FR	
2015-0308	8 FR	
2015-0310	10 FR	
2015-0312	12 FR	- Middle perforation
2015-0314	14 FR	- 1100 mm
2015-0316	16 FR	- X-ray opaque line
2015-0318	18 FR	
2015-0320	20 FR	




2015-02XX Silicone Round Drain




2015-03XX Silicone Round Drain

Round Drain Kit:

REF. No.	Size	Description
2015-9208	8 FR	8 FR end perf. Round, attached Trocar
2015-9210	10 FR	10 FR end perf. Round, attached Trocar
2015-9214	14 FR	14 FR end perf. Round, attached Trocar
2015-9216	16 FR	16 FR end perf. Round, attached Trocar
2015-9218	18 FR	18 FR end perf. Round, attached Trocar
2015-9220	20 FR	20 FR end perf. Round, attached Trocar
2015-9308	8 FR	8 FR mid. perf. Round, attached Trocar
2015-9310	10 FR	10 FR mid. perf. Round, attached Trocar
2015-9314	14 FR	14 FR mid. perf. Round, attached Trocar
2015-9316	16 FR	16 FR mid. perf. Round, attached Trocar
2015-9318	18 FR	18 FR mid. perf. Round, attached Trocar
2015-9320	20 FR	20 FR mid. perf. Round, attached Trocar



2015-92XX
Silicone Round Drain Kit




2015-93XX
Silicone Round Drain Kit

SILICONE ROUND DRAIN

Round Drain Set*:


REF. No.	Size	Description
0610-0208	8 FR	8 FR end perf. Round, 100 cc Reservoir
0610-0210	10 FR	10 FR end perf. Round, 100 cc Reservoir
0610-0212	12 FR	12 FR end perf. Round, 100 cc Reservoir
0610-0214	14 FR	14 FR end perf. Round, 100 cc Reservoir
0610-0216	16 FR	16 FR end perf. Round, 100 cc Reservoir
0610-0218	18 FR	18 FR end perf. Round, 100 cc Reservoir
0610-0220	20 FR	20 FR end perf. Round, 100 cc Reservoir
0610-0308	8 FR	8 FR mid. perf. Round, 100 cc Reservoir
0610-0310	10 FR	10 FR mid. perf. Round, 100 cc Reservoir
0610-0312	12 FR	12 FR mid. perf. Round, 100 cc Reservoir
0610-0314	14 FR	14 FR mid. perf. Round, 100 cc Reservoir
0610-0316	16 FR	16 FR mid. perf. Round, 100 cc Reservoir
0610-0318	18 FR	18 FR mid. perf. Round, 100 cc Reservoir
0610-0320	20 FR	20 FR mid. perf. Round, 100 cc Reservoir




0610-02XX Silicone Round Drain Set

Round Drain System*:


REF. No.	Size	Description
0610-9208	8 FR	8 FR end perf. Round, Trocar, 100 cc Reservoir
0610-9210	10 FR	10 FR end perf. Round, Trocar, 100 cc Reservoir
0610-9214	14 FR	14 FR end perf. Round, Trocar, 100 cc Reservoir
0610-9216	16 FR	16 FR end perf. Round, Trocar, 100 cc Reservoir
0610-9218	18 FR	18 FR end perf. Round, Trocar, 100 cc Reservoir
0610-9220	20 FR	20 FR end perf. Round, Trocar, 100 cc Reservoir
0610-9308	8 FR	8 FR mid. perf. Round, Trocar, 100 cc Reservoir
0610-9310	10 FR	10 FR mid. perf. Round, Trocar, 100 cc Reservoir
0610-9314	14 FR	14 FR mid. perf. Round, Trocar, 100 cc Reservoir
0610-9316	16 FR	16 FR mid. perf. Round, Trocar, 100 cc Reservoir
0610-9318	18 FR	18 FR mid. perf. Round, Trocar, 100 cc Reservoir
0610-9320	20 FR	20 FR mid. perf. Round, Trocar, 100 cc Reservoir



0610-92XX Silicone Round Drain System



0610-03XX Silicone Round Drain Set




0610-93XX Silicone Round Drain System

SILICONE ROUND DRAIN


Round Drain Bag System*:

REF. No.	Size	Description
0910-920820	8 FR	8 FR end perf. Round, Trocar, 100 cc Reservoir, 200 cc Bag
0910-921020	10 FR	10 FR end perf. Round, Trocar, 100 cc Reservoir, 200 cc Bag
0910-921420	14 FR	14 FR end perf. Round, Trocar, 100 cc Reservoir, 200 cc Bag
0910-921620	16 FR	16 FR end perf. Round, Trocar, 100 cc Reservoir, 200 cc Bag
0910-921820	18 FR	18 FR end perf. Round, Trocar, 100 cc Reservoir, 200 cc Bag
0910-922020	20 FR	20 FR end perf. Round, Trocar, 100 cc Reservoir, 200 cc Bag
0910-930820	8 FR	8 FR mid. perf. Round, Trocar, 100 cc Reservoir, 200 cc Bag
0910-931020	10 FR	10 FR mid. perf. Round, Trocar, 100 cc Reservoir, 200 cc Bag
0910-931420	14 FR	14 FR mid. perf. Round, Trocar, 100 cc Reservoir, 200 cc Bag
0910-931620	16 FR	16 FR mid. perf. Round, Trocar, 100 cc Reservoir, 200 cc Bag
0910-931820	18 FR	18 FR mid. perf. Round, Trocar, 100 cc Reservoir, 200 cc Bag
0910-932020	20 FR	20 FR mid. perf. Round, Trocar, 100 cc Reservoir, 200 cc Bag

* Round Drain Set, Round Drain System and Round Drain Bag System are also available with 150 cc, 200 cc or 400 cc Reservoir. Round Drain Bag System is also available with 400 cc or 800 cc Collection Bag.



0910-92XX20 Silicone Round Drain Bag System





0910-93XX20 Silicone Round Drain Bag System

SILICONE VAC DRAIN

Fortune® Silicone Vac Drain Features:

- Medical grade silicone for superior biocompatibility.
- X-ray opaque effectiveness.
- Multiple narrow ducts structure instead of perforation.
- Reduce withdrawal difficulty against tissue ingrowth.
- Rigid structure offers excellent kink-resistant character.
- Clear connecting tube provides inspection and drainage visibility.
- Available with 6 sizes: 4 round and 2 flat types.
- Sterilized double package.
- Available with pre-connected Trocar: Vac Drain Kit.

- A** Vac Drain
- B** Vac Drain Kit
- C** Vac Drain Set*
- D** Vac Drain System*
- E** Vac Drain Bag System*




SILICONE VAC DRAIN


SPECIFICATIONS

Vac Drain:

REF. No.	Size	Description
2014-0410	10 FR	- Full round Vac Drain
2014-0415	15 FR	- 200 mm round duct drain, 600 mm round tube
2014-0419	19 FR	- X-ray opaque effective
2014-0421	21 FR	- Full flat Vac Drain
2014-0408	3.5 x 8 mm	- 200 mm flat duct drain, 600 mm round tube
2014-0411	4 x 11 mm	- X-ray opaque effective




2014-0410/15/19/21




2014-0408/11

Vac Drain Kit:

REF. No.	Size	Description
2014-9410	10 FR	10 FR full Vac, attached Trocar
2014-9415	15 FR	15 FR full Vac, attached Trocar
2014-9419	19 FR	19 FR full Vac, attached Trocar
2014-9421	21 FR	21 FR full Vac, attached Trocar
2014-9408	3.5 x 8 mm	8 mm full Vac, attached Trocar
2014-9411	4 x 11 mm	11 mm full Vac, attached Trocar



2014-9410/15/19/21




2014-9408/11


SILICONE VAC DRAIN

Vac Drain Set*:

REF. No.	Size	Description
0610-0410	10 FR	10 FR full Vac, 100 cc Reservoir
0610-0415	15 FR	15 FR full Vac, 100 cc Reservoir
0610-0419	19 FR	19 FR full Vac, 100 cc Reservoir
0610-0421	21 FR	21 FR full Vac, 100 cc Reservoir
0610-0408	3.5 x 8 mm	8 mm full Vac, 100 cc Reservoir
0610-0411	4 x 11 mm	11 mm full Vac, 100 cc Reservoir




0610-0410/15/19/21




0610-0408/11

Vac Drain System*:

REF. No.	Size	Description
0610-9410	10 FR	10 FR full Vac, Trocar, 100 cc Reservoir
0610-9415	15 FR	15 FR full Vac, Trocar, 100 cc Reservoir
0610-9419	19 FR	19 FR full Vac, Trocar, 100 cc Reservoir
0610-9421	21 FR	21 FR full Vac, Trocar, 100 cc Reservoir
0610-9408	3.5 x 8 mm	8 mm full Vac, Trocar, 100 cc Reservoir
0610-9411	4 x 11 mm	11 mm full Vac, Trocar, 100 cc Reservoir



0610-9410/15/19/21






0610-9408/11

SILICONE VAC DRAIN


Vac Drain Bag System *:

REF. No.	Size	Description
0910-941020	10 FR	10 FR full Vac, Trocar, 100 cc Reservoir, 200 cc Bag
0910-941520	15 FR	15 FR full Vac, Trocar, 100 cc Reservoir, 200 cc Bag
0910-941920	19 FR	19 FR full Vac, Trocar, 100 cc Reservoir, 200 cc Bag
0910-942120	21 FR	21 FR full Vac, Trocar, 100 cc Reservoir, 200 cc Bag
0910-940820	3.5 x 8 mm	8 mm full Vac, Trocar, 100 cc Reservoir, 200 cc Bag
0910-941120	4 x 11 mm	11 mm full Vac, Trocar, 100 cc Reservoir, 200 cc Bag

* Vac Drain Set, Vac Drain System and Vac Drain Bag System are also available with 150 cc, 200 cc or 400 cc Reservoir. Vac Drain Bag System is also available with 400 cc or 800 cc Collection Bag.



0910-94(10/15/19/21)20



0910-94(08/11)20

SILICONE ENDO DRAIN

Fortune® Silicone Endo Drain Features:

- Medical grade silicone for superior biocompatibility.
- X-ray opaque line.
- Multiple narrow ducts structure instead of perforation.
- Reduce withdrawal difficulty against tissue ingrowth.
- Rigid structure offers excellent kink-resistant character.
- Clear connecting tube provides inspection and drainage visibility.
- Available with 4 sizes.
- Sterilized double package.
- Available with pre-connected Trocar: Endo Drain Kit.


- A Endo Drain**
- B Endo Drain Kit**
- C Endo Drain Set***
- D Endo Drain System***
- E Endo Drain Bag System***



SPECIFICATIONS

Endo Drain:


REF. No.	Size	Description
2015-0610	10 FR	
2015-0615	15 FR	
2015-0619	19 FR	- Transparent - 1200 mm
2015-0621	21 FR	- X-ray opaque line
2015-0624	24 FR	
2015-0629	29 FR	



2015-06XX Silicone Endo Drain

Endo Drain Kit:

REF. No.	Size	Description
2015-9610	10 FR	10 FR Endo, attached Trocar
2015-9615	15 FR	15 FR Endo, attached Trocar
2015-9619	19 FR	19 FR Endo, attached Trocar
2015-9621	21 FR	21 FR Endo, attached Trocar




2015-96XX Silicone Endo Drain Kit

SILICONE ENDO DRAIN


Endo Drain Set*:

REF. No.	Size	Description
0610-0610	10 FR	10 FR Endo, 100 cc Reservoir
0610-0615	15 FR	15 FR Endo, 100 cc Reservoir
0610-0619	19 FR	19 FR Endo, 100 cc Reservoir
0610-0621	21 FR	21 FR Endo, 100 cc Reservoir



Endo Drain System *:


REF. No.	Size	Description
0610-9610	10 FR	10 FR Endo, Trocar, 100 cc Reservoir
0610-9615	15 FR	15 FR Endo, Trocar, 100 cc Reservoir
0610-9619	19 FR	19 FR Endo, Trocar, 100 cc Reservoir
0610-9621	21 FR	21 FR Endo, Trocar, 100 cc Reservoir



Drain Bag System*:

REF. No.	Size	Description
0910-961020	10 FR	10 FR Endo, Trocar, 100 cc Reservoir, 200 cc Bag
0910-961520	15 FR	15 FR Endo, Trocar, 100 cc Reservoir, 200 cc Bag
0910-961920	19 FR	19 FR Endo, Trocar, 100 cc Reservoir, 200 cc Bag
0910-962120	21 FR	21 FR Endo, Trocar, 100 cc Reservoir, 200 cc Bag

* Endo Drain Set, Endo Drain System and Endo Drain Bag System are also available with 150 cc, 200 cc or 400 cc Reservoir. Endo Drain Bag System is also available with 400 cc or 800 cc Collection Bag.



SILICONE ENDO DRAIN (BENDABLE TROCAR + J CONNECTOR)

Features:

J Connector:


J Connector is designed to connect Endo Drain with the size of 10FR/15FR/19FR/21FR and Innerspring Reservoir. It is made of PVC which is a Non-Toxic material, with a little bendable character to reduce the risk of fracture. The side of J Connector designed to connect with Endo Drain is tapering, making assembly and disassembly more easily.



SPECIFICATIONS

Endo drain + J Connector:


2015-0610J	10FR	10 FR Endo, attached J Connector
2015-0615J	15FR	15 FR Endo, attached J Connector
2015-0619J	19FR	19 FR Endo, attached J Connector
2015-0621J	21FR	21 FR Endo, attached J Connector



2015-06xxJ

Endo Drain Kit (Trocar + J Connector):

2015-9610J	10FR	10 FR Endo, attached Trocar, J Connector
2015-9615J	15FR	15 FR Endo, attached Trocar, J Connector
2015-9619J	19FR	19 FR Endo, attached Trocar, J Connector
2015-9621J	21FR	21 FR Endo, attached Trocar, J Connector




2015-96xxJ

SILICONE ENDO DRAIN (BENDABLE TROCAR + J CONNECTOR)

Endo Drain with Bendable Trocar:


2015-9610B	10FR	10 FR Endo, attached Bendable Trocar
2015-9615B	15FR	15 FR Endo, attached Bendable Trocar
2015-9619B	19FR	19 FR Endo, attached Bendable Trocar
2015-9621B	21FR	21 FR Endo, attached Bendable Trocar



2015-96xxB

Endo Drain Kit (Bendable Trocar + J Connector):

2015-9610JB	10FR	10 FR Endo, attached Bendable Trocar, J Connector
2015-9615JB	15FR	15 FR Endo, attached Bendable Trocar, J Connector
2015-9619JB	19FR	19 FR Endo, attached Bendable Trocar, J Connector
2015-9621JB	21FR	21 FR Endo, attached Bendable Trocar, J Connector




2015-96xxJB

Bendable Trocar:

- Made of stainless steel.
- Bendable, sharp, smooth and burr free.
- Protection cover.
- Customized specification is available upon request.

2015-9110	10FR	
2015-9116	16FR	- Bendable
		- Stainless steel
2015-9119	19FR	- Protection cover
2015-9121	21FR	



2015-91xx

SILICONE HUBLESS ENDO DRAIN

Fortune® Silicone Hubless Endo Drain Features:

- Medical grade silicone for superior biocompatibility.
- Hubless design
- X-ray opaque line.
- Multiple narrow ducts structure instead of perforation.
- Reduce withdrawal difficulty against tissue ingrowth.
- Rigid structure offers excellent kink-resistant character.
- Clear connecting tube provides inspection and drainage visibility.
- Available with 4 sizes.
- Sterilized double package.
- Available with pre-connected Trocar: Hubless Endo Drain Kit.


- A** Hubless Endo Drain
- B** Hubless Endo Drain Kit
- C** Hubless Endo Drain Set*
- D** Hubless Endo Drain System*
- E** Hubless Endo Drain Bag System*



SPECIFICATIONS

Hubless Endo Drain:


REF. No.	Size	Description
2015-0710	10 FR	
2015-0715	15 FR	
2015-0719	19 FR	- Transparent - 1200 mm
2015-0721	21 FR	- X-ray opaque line
2015-0724	24 FR	
2015-0729	29 FR	



2015-07XX

Hubless Endo Drain Kit:

REF. No.	Size	Description
2015-9710	10 FR	10 FR Hubless Endo, attached Trocar
2015-9715	15 FR	15 FR Hubless Endo, attached Trocar
2015-9719	19 FR	19 FR Hubless Endo, attached Trocar
2015-9721	21 FR	21 FR Hubless Endo, attached Trocar




2015-97XX

SILICONE HUBLESS ENDO DRAIN

Hubless Endo Drain Set*:


REF. No.	Size	Description
0610-0710	10 FR	10 FR Hubless Endo, 100 cc Reservoir
0610-0715	15 FR	15 FR Hubless Endo, 100 cc Reservoir
0610-0719	19 FR	19 FR Hubless Endo, 100 cc Reservoir
0610-0721	21 FR	21 FR Hubless Endo, 100 cc Reservoir



0610-07XX

Hubless Endo Drain System*:

REF. No.	Size	Description
0610-9710	10 FR	10 FR Hubless Endo, Trocar, 100 cc Reservoir
0610-9715	15 FR	15 FR Hubless Endo, Trocar, 100 cc Reservoir
0610-9719	19 FR	19 FR Hubless Endo, Trocar, 100 cc Reservoir
0610-9721	21 FR	21 FR Hubless Endo, Trocar, 100 cc Reservoir




0610-97XX

Hubless Endo Drain Bag System*:

REF. No.	Size	Description
0910-971020	10 FR	10 FR Hubless Endo, Trocar, 100 cc Reservoir, 200 cc Bag
0910-971520	15 FR	15 FR Hubless Endo, Trocar, 100 cc Reservoir, 200 cc Bag
0910-971920	19 FR	19 FR Hubless Endo, Trocar, 100 cc Reservoir, 200 cc Bag
0910-972120	21 FR	21 FR Hubless Endo, Trocar, 100 cc Reservoir, 200 cc Bag

* Hubless Endo Drain Set, Hubless Endo Drain System and Hubless Endo Drain Bag System are also available with 150 cc, 200 cc or 400 cc Reservoir. Hubless Endo Drain Bag System is also available with 400 cc or 800 cc Collection Bag.



0910-97XX20

SILICONE HUBLESS ENDO DRAIN (BENDABLE TROCAR + J CONNECTOR)

Features:

J Connector:


J Connector is designed to connect Hubless Endo Drain with the size of 10FR/15FR/19FR/21FR and Innerspring Reservoir. It is made of PVC which is a Non-Toxic material, with a little bendable character to reduce the risk of fracture. The side of J Connector designed to connect with Hubless Endo Drain is tapering, making assembly and disassembly more easily.



SPECIFICATIONS

Hub-less Endo drain + J Connector:


2015-0710J	10FR	10 FR Endo, attached J Connector
2015-0715J	15FR	15 FR Endo, attached J Connector
2015-0719J	19FR	19 FR Endo, attached J Connector
2015-0721J	21FR	21 FR Endo, attached J Connector



2015-07xxJ

Hub-less Endo Drain Kit (Trocar + J Connector):

2015-9710J	10FR	10 FR Endo, attached Trocar, J Connector
2015-9715J	15FR	15 FR Endo, attached Trocar, J Connector
2015-9719J	19FR	19 FR Endo, attached Trocar, J Connector
2015-9721J	21FR	21 FR Endo, attached Trocar, J Connector




2015-97xxJ

SILICONE HUBLESS ENDO DRAIN (BENDABLE TROCAR + J CONNECTOR)

Hub-less Endo Drain with Bendable Trocar:


2015-9710B	10FR	10 FR Endo, attached Bendable Trocar
2015-9715B	15FR	15 FR Endo, attached Bendable Trocar
2015-9719B	19FR	19 FR Endo, attached Bendable Trocar
2015-9721B	21FR	21 FR Endo, attached Bendable Trocar



2015-97xxB

Hub-less Endo Drain Kit (Bendable Trocar + J Connector):

2015-9710JB	10FR	10 FR Endo, attached Bendable Trocar, J Connector
2015-9715JB	15FR	15 FR Endo, attached Bendable Trocar, J Connector
2015-9719JB	19FR	19 FR Endo, attached Bendable Trocar, J Connector
2015-9721JB	21FR	21 FR Endo, attached Bendable Trocar, J Connector



2015-97xxJB

Bendable Trocar:

- Made of stainless steel.
- Bendable, sharp, smooth and burr free.
- Protection cover.
- Customized specification is available upon request.

2015-9110	10FR	
2015-9116	16FR	- Bendable
		- Stainless steel
2015-9119	19FR	- Protection cover
2015-9121	21FR	



2015-91xx


ENDO CONNECTOR

Fortune® Endo Connector:




SPECIFICATIONS

REF. No.	Description
2015-8601	1:1 with plugs, connects 1 to 1 Endo-Drain of 19, 21, 24 & 29 FR, source port fits 10-13 mm tubing
2015-8603	3:1 with plugs, connects 1 to 3 Endo-Drain of 19, 21, 24 & 29 FR, source port fits 10-13 mm tubing



2015-8601



2015-8603



6 F., No. 29, Sec. 2, Jhongjheng E. Rd.,
Danshuei Dist, New Taipei City. 251 Taiwan
Tel.:(886) 2-2624-2233 Fax.:(886) 2-2624-2266
www.fortunemed.com Mail:info@fortunemed.com