

EU-MDR IEC ISO ANVISA

BalloonPass



Endopass



Quickpass



Hassonpass



Easypass



DISPOSABLE LAPAROSCOPIC TROCAR

Offering a variety of trocar products: complete specifications to meet diverse combination needs, excellent airtightness, and a comfortable feel.

- Comfortable Grip
- Complete Specifications
- Fully Detachable
- Stable and Reliable
- Double Sealing
- Excellent Airtightness
- Extra-Large Flow Channel
- Secure Fit

BalloonPass

Disposable Laparoscopic Trocar - Balloonpass




STABLE AND SECURE ATTACHMENT

Dual-layer fixation prevents trocar cannula slippage during surgery

BalloonPass

Secure fixation prevents slippage, while the smooth cannula reduces damage, bleeding, and air leakage. Dual sealing ensures cleanliness, and the detachable structure improves flexibility. Shallower insertion enhances the surgical field, and the large flow channel speeds pneumoperitoneum creation for safer, more efficient operations.




Blade-free Visual Tip High Transparency

High-transparency conical tip for clear view, gentle fiber separation, and accelerated healing.



Fixation Ring

Prevents trocar displacement, ensuring safety and sealing



Dynamic Sealing

Low resistance with excellent sealing for smooth operation




Durable Material

Scratch-resistant, tear-proof, and flexible for lasting performance





Rapid Airflow Design





Establishes pneumoperitoneum quickly, stabilizing pressure



Dual Sealing

Enhanced sealing with a cloverleaf valve to prevent leaks and ensure safety.

	Diameter: 5mm	Length: 100mm	TB 5×100
	Diameter: 10mm	Length: 100mm	TB 10x100
	Diameter: 11mm	Length: 100mm	TB 11x100
	Diameter: 12mm	Length: 100mm	TB 12x100
	Diameter: 15mm	Length: 100mm	TB 15x100

TCB 5X100 ϕ5mm		TCB 10X100 ϕ10mm	
TCB 11X100 ϕ11mm		TCB 12X100 ϕ12mm	

Endopass

Disposable Laparoscopic Trocar - Endopass




Endopass

Up to seven specifications with a transparent design for camera observation. Tapered tip reveals tissue layers, anti-fog system ensures clarity, and integrates with hospital camera systems. Reduces blind puncture risk, boosts success rate, lowers complications, suits obese and post-surgery patients, enhancing safety and confidence.



Blade-free Visual Tip High Transparency

High-transparency conical tip for clear view, gentle fiber separation, and accelerated healing.



Elastic Fins

Elastic fins with silicone wings close during insertion and open for secure, low-damage fixation



Dynamic Sealing

Low resistance with excellent sealing for smooth operation




Durable Material

Scratch-resistant, tear-proof, and flexible for lasting performance



Rapid Airflow Design


Establishes pneumoperitoneum quickly, stabilizing pressure



Dual Sealing

Enhanced sealing with a cloverleaf valve to prevent leaks and ensure safety.

	Diameter: 5mm	Length:75mm	TV 5x75
		Length:100mm	TV 5x100
		Length:150mm	TV 5x150
	Diameter: 10mm	Length:100mm	TV 10x100
		Length:150mm	TV 10x150
	Diameter: 11mm	Length:100mm	TV 11x100
		Length:150mm	TV 11x150
	Diameter: 12mm	Length:100mm	TV 12x100
		Length:150mm	TV 12x150

TCU 5X75	TCU 5X100	TCU 10X100	TCU 11X100	TCU 12X100
ϕ5mm	ϕ5mm	ϕ10mm	ϕ11mm	ϕ12mm
				
75mm	100mm	100mm	100mm	100mm


Easypass

Disposable Laparoscopic Trocar - Easypass




Easypass

Up to five specifications with a lightweight, user-friendly design. Sharp tip ensures smooth wall penetration, standard cannula provides stability, and sealing system maintains pneumoperitoneum. Cost-effective, compatible, easy to use, reduces prep time, and minimizes fatigue–ideal for laparoscopic surgeries.



Dual “V” Bladeless Tip

Smooth and sharp tapered transition enables easy fiber separation, reducing trauma and accelerating healing



Elastic Fins

Silicone wings close during insertion, open for secure, low-damage fixation



Dynamic Sealing

Low resistance with excellent sealing for smooth operation




Durable Material

Scratch-resistant, tear-proof, and flexible for lasting performance




Rapid Airflow Design

Establishes pneumoperitoneum quickly, stabilizing pressure



Dual Sealing

Enhanced sealing with a cloverleaf valve to prevent leaks and ensure safety.

	Diameter: 3mm	Length: 65mm	TS 3x65
	Diameter: 5mm	Length: 100mm	TS 5x100
	Diameter: 10mm	Length: 100mm	TS 10x100
	Diameter: 11mm	Length: 100mm	TS 11x100
	Diameter: 12mm	Length: 100mm	TS 12x100
	Diameter: 15mm	Length: 100mm	TS 15x100

Hassonpass

Disposable Laparoscopic Trocar - Hassonpass



HassonPass

Blunt tip prevents organ damage, tapered cannula ensures sealing, and fixation wings provide support. Suture anchors and wire posts secure the cannula, ideal for patients with prior abdominal surgeries.



Blunt-Tip Design

Reduces accidental puncture risk and organ or vascular injury by 75%, ideal for patients with prior abdominal surgeries.



Adjustable Tapered Fixation Wing System

3djustable tapered wings provide 360° support, adapt to wall thickness, ensure sealing, and stabilize pneumoperitoneum pressure.



Dynamic Sealing

Low resistance with excellent sealing for smooth operation



Durable Material

Scratch-resistant, tear-proof, and flexible for lasting performance



Rapid Airflow Design

Establishes pneumoperitoneum quickly, stabilizing pressure



Dual Sealing

Enhanced sealing with a cloverleaf valve to prevent leaks and ensure safety

	Diameter: 10mm	Comprimento: 100mm	TH 10×100
	Diameter: 11mm	Comprimento: 100mm	TH 11x100
	Diameter: 12mm	Comprimento: 100mm	TH 12x100


Quickpass

Disposable Laparoscopic Trocar - Quickpass




Quickpass

Up to four specifications with precision cutting, sharp triangular blades, and retractable protection. Efficient puncture, clean incisions, reduced strain, and streamlined design enhance safety and speed.



Automatic Puncture Protection System

Equipped with a protective sheath prevents blade damage and accidental organ injury



Elastic Fins

Silicone wings close during insertion, open for secure, low-damage fixation



Dynamic Sealing

Low resistance with excellent sealing for smooth operation



Durable Material

Scratch-resistant, tear-proof, and flexible for lasting performance



Rapid Airflow Design

Establishes pneumoperitoneum quickly, stabilizing pressure



Dual Sealing

Enhanced sealing with a cloverleaf valve to prevent leaks and ensure safety

	Diameter: 5mm	Length: 100mm	TP 5x100
	Diameter: 10mm	Length:100mm	TP 10x100
	Diameter: 11mm	Length: 100mm	TP 11x100
	Diameter: 12mm	Length: 100mm	TP 12x100